INSTRUCTIONS FOR USE

**Suture Aid Booties
Non-Sterile**

Aspen Surgical Products Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316
Phone 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

EC REP
Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands
Made in USA

English

Intended use:

Suture Aid Booties are intended to protect , tag and help locate sutures

Intended User/ Patient Target Groups:

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Instructions for use

1. After removing backing from cartridge, apply to any appropriate surface.
2. Insert tip of appropriate instrument into Suture Aid Bootie and remove from cartridge.

Extra Suture Aid Bootie

- Apply Extra Aid Booties over tip of appropriate instrument.

Precaution:

Device is single use only. Reuse of device could result in infection/ contamination and/ or device failure which could lead to patient harm. Device is intended for limited contact duration, less than or equal to 8 hours.

Sterilization Instructions

1. Place device into a kit or sterilization tray.
2. Using Ethylene Oxide sterilization method, sterilize according to the following parameters:

Cycle Parameter	Specification	Duration
Preconditioning		
Temperature	110°F +/- 10°F	
Humidity	65% +/- 15%	10.35 Hrs - 72 Hrs
Chamber		
Vacuum	21" +/- 1.0 inHgA	
Humidification/Dwell	22.5" +/- 1.0 HgA	60 +30/-0 minutes
Gas Inject/Dwell	10.3" +/- 0.5 Hg rise	360 +5/-0 minutes
Jacket Temperature	135°F	
Aeration		
Temperature	120°F	24 Hrs - 72 Hrs Max

Disposal:

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Notice to Users and/or Patients in EU:

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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