INSTRUCTIONS FOR USE

Disposable Light Handle Covers – Non-Sterile

MD R_x ONLY

Aspen Surgical Products Inc. 6945 Southbelt Drive SE Caledonia, MI 49316 Phone 616-698-7100 Toll-Free 888-364-7004 aspensurgical.com	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands	Flexible Light Handle Covers-Made in China
		Rigid Light Handle Covers-Made in USA

English

Intended use:

intended to provide a sterile interface between the operating room surgical lamps (non-sterile) and the Surgeon or Nurse (sterile).

Intended User/ Patient Target Groups:

Device is non-patient contacting and is to be used by Healthcare Professionals.

Provides a means of manipulating O.R. Lights.

Precautions:

- Verify securement of the light handle cover to the adaptor.
- Do not use if product appears visually damaged or degraded.
- This device is single use only. Reuse of device could result in infection/ contamination and/ or device failure which could lead to patient harm.

Sterilization Instructions

- Place device into a kit or sterilization tray.
- Using Ethylene Oxide sterilization method, sterilize according to the following parameters:

Cycle Parameter	Specification	Duration
Preconditioning		
Temperature	110°F +/- 10°F	
Humidity	65% +/- 15%	10.35 Hrs - 72 Hrs
Chamber		
Vacuum	21" +/- 1.0 inHgA	
Humidification/Dwell	22.5" +/- 1.0 HgA	60 +30/-0 minutes
Gas Inject/Dwell	10.3" +/- 0.5 Hg rise	360 +5/-0 minutes
Jacket Temperature	135°F	
Aeration		
Temperature	120°F	24 Hrs - 72 Hrs Max

Dispose:

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Notice to Users and/or Patients in EU:

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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