## INSTRUCTIONS FOR USE

# O.R. Camera Cover

Non-Sterile

MD R<sub>x</sub> ONLY

EC REP

 Emergo Europe  
 Westervoortseidijk 60  
 6827 AT Arnhem  
 The Netherlands

For Country of Origin See Product Label

 Aspen Surgical Products Inc.  
 6945 Southbelt Drive SE  
 Caledonia, MI 49316  
 Phone +1 616-698-7100  
 Toll-Free 888-364-7004  
 aspensurgical.com

## English

## Intended Use

Intended to provide a sterile interface between the operating room camera (non-sterile) and the Surgeon or Nurse (sterile).

## Intended User/ Patient Target Groups

Device is non-patient contacting and is to be used by Healthcare Professionals. Provides a means of manipulating O.R. Camera.

## Instruction for Use:

- Verify securement of the camera cover to the adaptor.
- Device is intended for use less than or equal to 16 hours.

## Sterilization Instructions:

- Place device into a kit or sterilization tray.
- Using Ethylene Oxide sterilization method, validate according to ISO 11135, sterilize according to the following parameters:

Cycle Parameter	Specification	Duration
<b>Preconditioning</b>		
Temperature	110°F +/- 10°F	
Humidity	65% +/- 15% RH	10.35 Hrs - 72 Hrs
<b>Chamber</b>		
Vacuum	21" +/- 1.0 inHgA	
Humidification/Dwell	22.5" +/- 1.0 HgA	60 +30/-0 minutes
Gas Inject/Dwell	10.3" +/- 0.5 Hg rise	360 +5/-0 minutes
Jacket Temperature	135° +/- 5° F	
<b>Aeration*</b>		
Temperature	120° +/- 10° F	24 Hrs - 72 Hrs Max

\*Proper aeration is essential to eliminating the risk of residuals to meet ISO 10993-7 requirements.

## Warning

Reuse of device could result in infection/contamination and/or device failure which could lead to patient harm.

## Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

## Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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### AB'deki Kullanıcılar ve/veya Hastalara Yönelik Bildirim
Emergo Consulting (UK) Limited