Vessel Loops

aspensurgical.

Non-Sterile

Aspen Surgical Products, Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316 USA
Phone +1 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

EMERGO
Emergo Europe
Westervoortseedijk 60
6822 AT Arnhem
The Netherlands

CE
2797

For Country of Origin See Product Label

English

Intended Use

Vessel loops are intended to be used to occlude, retract, and identify arteries, veins, tendons, and nerves in surgical procedures.

Intended User/Patient Target Groups

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Instruction for Use

Device is intended for limited contact duration, less than or equal to 8 hours.

Sterilization Instructions:

- Place device into a kit or sterilization tray.
- Using Ethylene Oxide sterilization method, validate according to ISO 11135, sterilize according to the following parameters:

Cycle Parameter	Specification	Duration
Preconditioning		
Temperature	110°F +/- 10°F	
Humidity	65% +/- 15% RH	10.35 Hrs - 72 Hrs
Chamber		
Vacuum	21" +/- 1.0 inHgA	
Humidification/Dwell	22.5" +/- 1.0 HgA	60 +30/-0 minutes
Gas Inject/Dwell	10.3" +/- 0.5 Hg rise	360 +5/-0 minutes
Jacket Temperature	135° +/- 5° F	
Aeration*		
Temperature	120° +/- 10° F	24 Hrs - 72 Hrs Max

*Proper aeration is essential to eliminating the risk of residuals to meet ISO 10993-7 requirements.

Warning

To help reduce the potential for infection and/or other complications, do not re-use

Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional and/or local laws and regulations.

Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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AB'deki Kullanıcılar ve/veya Hastalara Yönelik Bildirim
Emergo Consulting (UK) Limited