

Bard-Parker® SafeSwitch® Disposable Scalpel Handle & Disposable Scalpel Handle Cover

Intended Use

Intended to be used with surgical blades for tissue separation and other procedures that require a sharp surgical blade to puncture or cut.

Key Features

- Helps provide protection from injuries that may occur before use, during steps of a procedure, or after use
- Designed for both left and right-handed surgeons
- Easy one-step blade removal design
- Features a retractable shield that boasts easy one-handed activation
- Not made with natural rubber latex
- Disposable

Scalpel Handle


- Compatible with the following Bard-Parker® blade sizes: 10, 11, 15 and 15C (blades not included)
- Designed with a similar weight to existing Bard-Parker® conventional and protected blade system products
- Compatible with both Bard-Parker® stainless steel and carbon Rib Back® blade styles

Scalpel Handle Cover

- Compatible with Bard-Parker® conventional size 3 handle (handle not included)

Bard-Parker® SafeSwitch® Disposable Scalpel Handle & Disposable Scalpel Handle Cover

Intended to be used with surgical blades for tissue separation and other procedures that require a sharp surgical blade to puncture or cut.



Blades sold separately.

Technical Specifications

Part No.	Description	Size	Sterility	Qty
ST-930D	Bard-Parker® SafeSwitch® Disposable Scalpel Handle	10, 11, 15, 15C	Sterile	8/Box, 32/Case
ST-830D	Bard-Parker® SafeSwitch® Disposable Scalpel Handle Cover	3	Sterile	8/Box, 32/Case

NEW!

Bard-Parker® and SafeSwitch® are registered trademarks of Aspen Surgical. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2017 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1019934 rev 2 10-FEB-2017 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

Phone 616-698-7100
Toll-free 888-364-7004

www.aspensurgical.com

Enhancing outcomes for patients and their caregivers:

Hill-Rom