



# SurgiSafe® SelectMat™ Perforated Absorbent Floor Mat

## Intended Use

Used to help keep floors in the healthcare setting clean, dry, and slip-free particularly during fluid intense procedures.

## Key Features

- Used to manage infectious waste fluid generated during wet procedures in surgical settings
- Enables quick O.R. turnover time by retaining the fluid generated in the procedure
- Quickly absorbs irrigants, blood, and bodily fluids through superior wicking, absorbency, and fluid retention properties
- Aids facilities in compliance with OSHA Bloodborne Pathogens and General Environmental Controls Standards<sup>1</sup>
- Perforations every 15 inches allows user to select the desired length
- Available with or without fluid barrier backing
- Provides a durable surface for staff to stand on
- Strong and dense, allowing staff to roll over them without catching the mat in their casters
- Non-Sterile, intended for single-use, and not made with natural rubber latex

# SurgiSafe® SelectMat™ Perforated Absorbent Floor Mat

Used to help keep floors in the healthcare setting clean, dry, and slip-free particularly during fluid intense procedures.



## Technical Specifications

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
81575P	Yellow	Medium	Yes	15" x 75'	1 Roll/Box
81575	Yellow	Medium	No	15" x 75'	1 Roll/Box
83075P	Yellow	Medium	Yes	30" x 75'	1 Roll/Box
83075	Yellow	Medium	No	30" x 75'	1 Roll/Box

<sup>1</sup> OSHA. Occupational Safety and Health Standards. Toxic and Hazardous Substances. Bloodborne Pathogens. 1910.1030.

SurgiSafe® is a registered trademark of Aspen Surgical Products, Inc. SelectMat™ is a trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2018 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1036223 rev 1 13-JUN-2018 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

**Phone 616-698-7100**  
**Toll-free 888-364-7004**

**[www.aspensurgical.com](http://www.aspensurgical.com)**

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**