



AquaNot® Floor Suction Device

Intended Use


Intended to help control fluids generated during wet procedures by maneuvering the AquaNot® Floor Suction Device with a push of your foot to the areas affected by fluids.

Key Features

- The AquaNot® Device has similar characteristics as the WaterBug Quiet Floor Suction® Device, but has a lower profile
- Efficiently suction fluid waste to help reduce risks and maintains a dry operating room floor
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- The raised design allows for easy repositioning, controlling fluids over a large floor area
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-Sterile, intended for single-use, and not made with natural rubber latex

AquaNot® Floor Suction Device

Intended to help control fluids generated during wet procedures by maneuvering the AquaNot® Floor Suction Device with a push of your foot to the areas affected by fluids.



Technical Specifications

Part No.	Description	Tubing	Qty
97020	AquaNot® Floor Suction Device	Yes	20/Case
97060	AquaNot® Floor Suction Device	Yes	60/Case
97050-NT	AquaNot® Floor Suction Device	No	50/Case

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2015 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1014393 rev 1 03-MAR-2015 ENG - US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

Phone 616-698-7100

Toll-free 888-364-7004

www.aspensurgical.com