



Oxboro[®] Corner Protectors

Intended Use

Intended to help protect instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays.

Key Features

- Aspen Surgical's corner protectors are designed to create a flexible barrier between the tray's sharp corners and the sterilization wrap
- They can be used with steam, gas, ozone, and low temperature plasma sterilization
- Reusable
- Not made with natural rubber latex

Oxboro® Corner Protectors

Intended to help protect instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays.



Technical Specifications

Part No.	Description	Sterility	Qty/Bag
CP-80NS	Corner Protectors	Non-Sterile	80

Oxboro® is a registered trademark of Aspen Surgical. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2015 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1000885 rev 1 16-APR-2015 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

Phone 616-698-7100

Toll-free 888-364-7004

www.aspensurgical.com