

## ViraSorb® Solidifier

- Evenly gels canister contents in seconds without stirring or agitating
- Non-toxic, glutaraldehyde-free, and can be safely incinerated
- Does not revert back to a fluid state, helping to eliminate spillage and leakage, qualifying it as Solid Regulated Medical Waste
- Pour nozzles are included with multi-use bottles
- Non-Sterile and not made with natural rubber latex



Part No.	Description	Solidifies	Qty
1000	ViraSorb® Solidifier – Single-Use Bottle	1,200cc	96/Case
1001	ViraSorb® Solidifier – Single-Use Bottle	1,500cc	96/Case
1002	ViraSorb® Solidifier – Single-Use Bottle	2,000cc	36/Case
1003	ViraSorb® Solidifier – Single-Use Bottle	3,000cc	36/Case
1005	ViraSorb® Solidifier – Single-Use Pouch	500cc	200/Case
1010	ViraSorb® Solidifier – Multi-Use Bottle	10,000cc	12/Case
1020	ViraSorb® Solidifier – Multi-Use Bottle	20,000cc	12/Case

AquaNot®, CinchMat®, ErgoSupport®, EvacMat®, SurgiSafe®, ViraSorb®, WaterBoom®, and WaterBug Quiet Floor Suction® are registered trademarks of Aspen Surgical Products Inc. Colby™ and WaterBug™ are trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2015 Aspen Surgical, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1006928 rev 4 16-MAR-2015 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

**Phone 616-698-7100**  
**Toll-free 888-364-7004**

**[www.aspensurgical.com](http://www.aspensurgical.com)**

## Aspen Surgical®

Floor Safety and Fluid Management



Enhancing outcomes for patients and their caregivers:

**Hill-Rom**

## Colby™ Fluid Management Products

Our Colby™ fluid management product line provides a convenient and effective means of managing infectious waste fluid generated during wet procedures. Our products help to reduce the likelihood of injuries associated with wet floors. They aid in protecting healthcare professionals and housekeeping personnel from contact with contaminated substances and help support your efforts to comply with OSHA's Bloodborne Pathogens Standard. Our products do not change the contents of medical waste, but rather help your team to more effectively contain it.

### SurgiSafe® Absorbent Floor Mats

- Helps keep floors clean, dry, and slip-free, even during fluid intense procedures
- Absorbs irrigants, blood, and bodily fluids through unique wicking, absorbency, and fluid retention properties, while conforming to uneven floor surfaces without bunching
- Helps protect your staff from slips, trips, and falls in the operating room
- Helps reduce staff exposure to infectious waste
- Helps eliminate re-using blankets and towels on patients
- Helps reduce O.R. turnaround time
- Contains fluid in one area for easy clean up
- Non-Sterile, intended for single-use, and not made with natural rubber latex

### SurgiSafe® Standard Absorbent Floor Mats

- Made with a minimum of 60% pre-consumer recycled material
- Available in three absorbency levels as well as with and without fluid barrier backing for customer convenience


Pink – High level of fluid absorbency, will absorb 16x's its dry weight  
 White – Medium level of fluid absorbency, will absorb 14x's its dry weight  
 Green – Low level of fluid absorbency, will absorb 11x's its dry weight

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
83430	Pink	High	No	28" x 40"	30/Case
83472	Pink	High	No	28" x 72"	30/Case
83330	Pink	High	Yes	28" x 40"	30/Case
83372	Pink	High	Yes	28" x 72"	30/Case
81860	White	Medium	No	18" x 40"	60/Case
818100	White	Medium	No	18" x 18"	100/Case
81850	White	Medium	No	18" x 18"	50/Case
83610	White	Medium	No	28" x 40"	10/Case
83630	White	Medium	No	28" x 40"	30/Case
83625-L	White	Medium	No	28" x 72"	12/Case
83672	White	Medium	No	28" x 72"	30/Case
83530	White	Medium	Yes	28" x 40"	30/Case
83572	White	Medium	Yes	28" x 72"	30/Case
84030	Green	Low	No	28" x 40"	30/Case
84072	Green	Low	No	28" x 72"	30/Case
840125	Green	Low	No	28" x 125"	24/Case
83030	Green	Low	Yes	28" x 40"	30/Case
83072	Green	Low	Yes	28" x 72"	30/Case
830125	Green	Low	Yes	28" x 125"	24/Case


### SurgiSafe® Specialty Absorbent Floor Mats

- Designed to include unique features beyond those offered by the Standard Absorbent Floor Mat Line
- Low profile
- Offered in large sizes
- Blue mat features fluid barrier backing to help keep floor beneath dry and clean
- Yellow mat is durable and absorbs fluid from the top and bottom
- CinchMat® – The blue mat with a draw cord for ease of disposal
- EvacMat® – The yellow mat with fluid suction capability

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
82350	Blue	Medium	Yes	23" x 40"	50/Case
84625	Blue	Medium	Yes	46" x 40"	25/Case
84610	Blue	Medium	Yes	46" x 40"	10/Case
84672	Blue	Medium	Yes	40" x 72"	10/Case
81860-2	Yellow	Medium	No	18" x 40"	60/Case
83630-2	Yellow	Medium	No	36" x 40"	30/Case
83610-2	Yellow	Medium	No	36" x 40"	10/Case
83615-2L	Yellow	Medium	No	36" x 72"	15/Case
83625-2L	Yellow	Medium	No	36" x 72"	12/Case
836125-2	Yellow	Medium	No	36" x 125"	10/Case
83610-M	White	Medium	No	36" x 40"	5/Case
	White	Medium	No	28" x 40"	5/Case
83630-M	Yellow	Medium	No	36" x 40"	15/Case
	White	Medium	No	28" x 40"	15/Case
84625-C	CinchMat®	Medium	Yes	46" x 40"	10/Case
95005	EvacMat®	Medium	Yes	36" x 48"	5/Case



SurgiSafe® Standard Mats



SurgiSafe® Specialty Mats



SurgiSafe® CinchMat®



SurgiSafe® EvacMat®



SurgiSafe® O.R. Cord Cover



ErgoSupport® Anti-Fatigue Mat

### SurgiSafe® O.R. Cord Cover

- Intended to help organize tubes and cords in the O.R.
- Serves as a "beacon" so healthcare personnel know where floor cords and tubes are located and helps decrease the chance of staff tripping on them
- Made of a durable material that can be rolled over
- Features a residue-free adhesive, helping it stay in one place
- Non-Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Color	Size	Qty
96050	O.R. Cord Cover	Yellow	9.5" x 24"	50/Case

### ErgoSupport® Anti-Fatigue Mat

- Intended to help reduce fatigue that healthcare personnel experience during long hours of standing in one location
- Uses a dual foam and gel technology to provide comfort under the foot
- The top surface provides a friction level with damp or dry surgical booties
- The bottom surface is made of a high-traction material to provide slip resistance
- Available in three sizes
- Water-resistant
- Reusable for up to 6-12 months
- Non-sterile and not made with natural rubber latex

Part No.	Description	Color	Size	Qty
81005	Anti-Fatigue Mat	Red	13" x 15"	5/Case
82001	Anti-Fatigue Mat	Red	18" x 24"	1/Case
82005	Anti-Fatigue Mat	Red	18" x 24"	5/Case
83002	Anti-Fatigue Mat	Red	20" x 48"	2/Case

## WaterBug Quiet Floor Suction® Device

- Efficiently suctions fluid waste to help reduce risks and maintain a dry operating room floor
- Its raised design and 12' of flexible tubing allows for easy repositioning
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-Sterile, intended for single-use, and not made with natural rubber latex

## AquaNot® Floor Suction Device

- Has similar characteristics as the WaterBug Quiet Floor Suction® Device, but has a lower profile
- Efficiently suctions fluid waste to help reduce risks and maintains a dry operating room floor
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Tubing	Qty
90010	WaterBug Quiet Floor Suction® Device	Yes	10/Case
90020	WaterBug Quiet Floor Suction® Device	Yes	20/Case
90030	WaterBug Quiet Floor Suction® Device	Yes	30/Case
97020	AquaNot® Floor Suction Device	Yes	20/Case
97060	AquaNot® Floor Suction Device	Yes	60/Case
97050-NT	AquaNot® Floor Suction Device	No	50/Case

## WaterBoom® Floor Suction Strip

- Creates a dam to help hold back, contain, and quietly suction an unlimited amount of waste fluid
- Provides strong suction, while conforming to your O.R. floor
- Its low profile helps you step over it without suction loss
- Its 12' tubing attaches from the side to help reduce tripping hazards
- This unobtrusive device does not require repositioning and is a good choice for those who prefer nothing at their feet
- Non-Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Length	Tubing	Qty
94010	WaterBoom® Floor Suction Strip	32"	Yes	10/Case
94030	WaterBoom® Floor Suction Strip	32"	Yes	30/Case