



Endoscopic Kittner

Intended Use

Intended for blunt dissection of tissue.

Key Features

- Our blunt dissecting instrument is 38cm long with a 5mm stainless steel shaft
- Available in both single and dual tip
- Radiopaque tip
- The brushed finish allows for safe use during laser surgery
- Sterile, intended for single-use, and not made with natural rubber latex

Endoscopic Kittner

Intended for blunt dissection of tissue.



Technical Specifications

| Part No. | Description | Sterility | Per/Box |
|----------|-------------|-----------|---------|
| KT-9101 | Single Tip | Sterile | 20 |
| KT-9105 | Dual Tip | Sterile | 20 |

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2015 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1014673 rev 1 14-APR-2015 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

Phone 616-698-7100

Toll-free 888-364-7004

www.aspensurgical.com