



# SurgiSafe® EvacMat® Absorbent Floor Mat

## Intended Use

Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor.

## Key Features

- Extremely quiet operation helps to decrease distraction of the O.R. staff – as quiet as an indoor conversation
- The top layer of the mat contains wicking to help prevent pooling
- Works well with wall suction
- Simple to use and clean up – red bag friendly
- Low profile allows for the O.R. staff to roll equipment over the mat without loss of suction
- Large enough for two people
- Provides cushion to help reduce fatigue during long procedures
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-Sterile, intended for single-use, and not made with natural rubber latex

## SurgiSafe® EvacMat® Absorbent Floor Mat

Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor.



### Technical Specifications

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
95005	Yellow	Medium	Yes	36" x 48"	5/Case

EvacMat® and SurgiSafe® are a registered trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2015 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1014391 rev 1 05-MAR-2015 ENG - US

For further information about this product or a service,  
please contact your local Aspen Surgical representative  
or visit our webpage:

**Phone 616-698-7100**  
**Toll-free 888-364-7004**

**[www.aspensurgical.com](http://www.aspensurgical.com)**