

# Specialty Disposable Light Handle

For use with Trumpf Medical™ ALC Plus Surgical Light  
Models: TruLight® 5000 and iLED® 3/5

## Intended Use

The disposable light handle is a single use sterile product. It is intended to be used as an accessory to surgical lamps. This disposable handle provides a sterile interface between the operating room surgical lamp (non-sterile) and the surgeon or scrub nurse (sterile).

## Key Features

- Designed to have a similar look and feel as the Trumpf Medical™ reusable handle
- Snap-on feature allows for quick and easy installation and removal
- Lens clarity is equivalent to the reusable handle, providing rapid response of the laser and sensor
- Light intensity can be adjusted with the touch sensors while the disposable light handle is being used
- Ultra-durable clear lens
- Avoids having to sterilize reusable handles
- Helps enable quick O.R. turnaround times
- Intended for single-use and not made with natural rubber latex

## Specialty Disposable Light Handle

For use with Trumpf Medical™ ALC Plus Surgical  
Light Models: TruLight® 5000 and iLED® 3/5



### Technical Specifications

Part No.	Description	Color	Sterility	Qty
LT-ALC01	Disposable Light Handle Cover	Green	Sterile	18/Box

Trumpf® and Trulight® are trademarks of Trumpf GmbH + Co. KG, Germany. iLED® is a registered trademark of TRUMPF Medizin Systeme GmbH + Co. KG. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2015 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1005864 rev 4 27-FEB-2015 ENG – US

For further information about this product or a service,  
please contact your local Aspen Surgical representative  
or visit our webpage:

**Phone 616-698-7100**

**Toll-free 888-364-7004**

**[www.aspensurgical.com](http://www.aspensurgical.com)**