



ViraSorb® Solidifier

Intended Use

Used to manage infectious waste fluid generated during wet procedures in surgical settings.

Key Features

- Evenly gels canister contents in seconds without stirring or agitating
- Non-toxic, glutaraldehyde-free, and can be safely incinerated
- It does not revert back to a fluid state, helping to eliminate spillage and leakage, qualifying it as a Solid Regulated Medical Waste
- Pour nozzles are included with multi-use-bottles
- Intended for single-use and not made with natural rubber latex

ViraSorb® Solidifier

Used to manage infectious waste fluid generated during wet procedures in surgical settings.



Technical Specifications

Part No.	Description	Solidifies	Qty/Case
1000	ViraSorb® Solidifier – Single-use Bottle	1,200cc	96
1001	ViraSorb® Solidifier – Single-use Bottle	1,500cc	96
1002	ViraSorb® Solidifier – Single-use Bottle	2,000cc	36
1003	ViraSorb® Solidifier – Single-use Bottle	3,000cc	36
1005	ViraSorb® Solidifier – Single-use Pouch	500cc	200
1010	ViraSorb® Solidifier – Multi-use Bottle	10,000cc	12
1020	ViraSorb® Solidifier – Multi-use Bottle	20,000cc	12

ViraSorb® is a registered trademark of Aspen Surgical. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2015 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1014679 rev 1 16-APR-2015 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

Phone 616-698-7100
Toll-free 888-364-7004

www.aspensurgical.com