



Standard Disposable Light Handle

Intended Use

Intended to provide a sterile interface between the operating room surgical light and the surgeon or nurse.

Please Note: An adaptor is needed to use the disposable light handles. Adaptors are available for most popular lights. Contact your local sales representative to begin the light handle conversion process today.

Key Features

- Two systems offered for convenience: Rigid or Flexible Covers
- Creates uniformity throughout hospital by utilizing one light handle product for multiple surgical lights
- Helps to save time and money associated with reprocessing reusable handles
- Convenient kit item, helps enable quicker O.R. turnaround times
- Single use
- Sterile and not made with natural rubber latex

Standard Disposable Light Handle

Creating a sterile interface in the O.R.



Flexible Handle



Rigid Handle

Technical Specifications

Part No.	Description	Color	Per/Box
LT-R01	Rigid Handle	Green	24
LT-F01B	Flexible Cover (1 each)	Green	100
LT-F02	Flexible Cover (2 each)	Green	80
LT-F03	Flexible Cover (3 each)	Green	80

Aspen Surgical reserves the right to make changes without notice in design, specifications and models. The only warranty Aspen Surgical makes is the express written warranty extended on the sale or rental of its products.

©2014 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-10008573 rev 1 23-JUL-2014 ENG - US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

Phone 616-698-7100

Toll-free 888-364-7004

www.aspensurgical.com