

Endoscopic Kittner

Intended Use


Intended for blunt dissection of tissue.

Key Features

- Our blunt dissecting instrument is 38cm long with a 5mm stainless steel shaft
- Available in both single and dual tip
- Radiopaque tip
- The brushed finish allows for safe use during laser surgery
- Sterile, intended for single-use, and not made with natural rubber latex

Endoscopic Kittner

Intended for blunt dissection of tissue.



Endoscopic Kittner Technical Specifications

Part No.	Description	Sterility	Qty
KT-9101	Single Tip	Sterile	20/Box
KT-9105	Dual Tip	Sterile	20/Box

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1014673 rev 3 8-Dec-2020 ENG - US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

