



SurgiSafe® SelectMat™ Perforated Absorbent Floor Mat

Intended Use

Used to help keep floors in the healthcare setting clean, dry, and slip-free particularly during fluid intense procedures.

Key Features

- Used to manage infectious waste fluid generated during wet procedures in surgical settings
- Enables quick O.R. turnover time by retaining the fluid generated in the procedure
- Quickly absorbs irrigants, blood, and bodily fluids through superior wicking, absorbency, and fluid retention properties
- Complies with OSHA's mandated exposure control plan for the worksite by reducing the likelihood of employee exposure to bloodborne pathogens ①
- Perforations every 15 inches allows user to select the desired length
- Available with or without fluid barrier backing
- Provides a durable surface for staff to stand on
- Strong and dense, allowing staff to roll over them without catching the mat in their casters
- Non-Sterile, intended for single-use, and not made with natural rubber latex

SurgiSafe® SelectMat™ Perforated Absorbent Floor Mat

Used to help keep floors in the healthcare setting clean, dry, and slip-free particularly during fluid intense procedures.



SurgiSafe® SelectMat™ Perforated Absorbent Floor Mat Technical Specifications

Part No.	Color	Absorbency	Fluid Barrier Level	Size Backing	Qty
81575P	Yellow	Medium	Yes	15" x 75'	1 Roll/Box
81575	Yellow	Medium	No	15" x 75'	1 Roll/Box
83075P	Yellow	Medium	Yes	30" x 75'	1 Roll/Box
83075	Yellow	Medium	No	30" x 75'	1 Roll/Box

SurgiSafe® is a registered trademark of Aspen Surgical Products, Inc. SelectMat™ is a trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

© OSHA, Occupational Safety and Health Standards. Toxic and Hazardous Substances, Bloodborne Pathogens, 1910.1030.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1036223 rev 2 9-Dec-2020 ENG - US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

