

# Sterion<sup>®</sup> Skin Closures

## Intended Use

Intended for skin closure and closure support for sutures and staples.

## Key Features

- Available in filament reinforced and flexible styles
- No adhesive enhancers required
- Less painful to apply than suturing or stapling, providing more comfort to the patient
- Sterile, intended for single-use and not made with natural rubber latex

## Filament Reinforced Skin Closures

- Reinforcing filaments offer added strength for large skin surface closures
- Filament Reinforced Skin Closure style is accepted as the standard by most healthcare professionals
- Breathable material designed to allow moisture transmission, to help reduce the chance of skin maceration
- Strips are curl-resistant

## Flexible Skin Closures

- Non-woven flexible backing helps reduce the chance of skin shear that can occur when edema or flexing puts excessive pressure on the skin closure ends

## Sterion® Skin Closures

Intended for skin closure and closure support for sutures and staples



### Technical Specifications – Filament

| Part No.  | Description                       | mm         | inches      | Strips/Pouch | Pouches/Box |
|-----------|-----------------------------------|------------|-------------|--------------|-------------|
| 751040pbx | Filament Reinforced Skin Closures | 3 x 75     | 1/8 x 3     | 5            | 50          |
| 751041pbx | Filament Reinforced Skin Closures | 6 x 75     | 1/4 x 3     | 3            | 50          |
| 751042pbx | Filament Reinforced Skin Closures | 6 x 38     | 1/4 x 1 1/2 | 6            | 50          |
| 751046pbx | Filament Reinforced Skin Closures | 6 x 100    | 1/4 x 4     | 10           | 50          |
| 751047pbx | Filament Reinforced Skin Closures | 12 x 100   | 1/2 x 4     | 6            | 50          |
| 751048pbx | Filament Reinforced Skin Closures | 25.4 x 127 | 1 x 5       | 4            | 25          |

### Technical Specifications – Flexible

| Part No.  | Description            | mm       | inches      | Strips/Pouch | Pouches/Box |
|-----------|------------------------|----------|-------------|--------------|-------------|
| 753040pbx | Flexible Skin Closures | 3 x 75   | 1/8 x 3     | 5            | 50          |
| 753041pbx | Flexible Skin Closures | 6 x 75   | 1/4 x 3     | 3            | 50          |
| 753042pbx | Flexible Skin Closures | 6 x 38   | 1/4 x 1 1/2 | 6            | 50          |
| 753046pbx | Flexible Skin Closures | 6 x 100  | 1/4 x 4     | 10           | 50          |
| 753047pbx | Flexible Skin Closures | 12 x 100 | 1/2 x 4     | 6            | 50          |
| 753048pbx | Flexible Skin Closures | 25 x 125 | 1 x 5       | 4            | 25          |

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2017 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1027754 rev 1 17-JAN-2017 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

**Phone 616-698-7100**  
**Toll-free 888-364-7004**

**[www.aspensurgical.com](http://www.aspensurgical.com)**

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**