



## SurgiSafe® CinchMat® Absorbent Floor Mat

### Intended Use


Combines the effective absorbency of the SurgiSafe® mat with an easier way to handle the mat during disposal.

### Key Features

- Pull cord design helps make clean-up and O.R. turnaround time simple and fast
- The top layer has wicking abilities that quickly absorbs irrigants, blood, and bodily fluids from the O.R. floor to help prevent pooling and splashing
- Skid-resistant fluid barrier backing helps ensure that the floor beneath the mat remains dry and prevents staining
- Large enough for two people
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-Sterile, intended for single-use, and not made with natural rubber latex

# SurgiSafe® CinchMat® Absorbent Floor Mat

Combines the effective absorbency of the SurgiSafe® mat with an easier way to handle the mat during disposal.



## SurgiSafe® CinchMat® Absorbent Floor Mat Technical Specifications

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
84625-C	Blue	Medium	Yes	46" x 40"	10/Case

CinchMat® and SurgiSafe® are registered trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1014392 rev 3 9-Dec-2020 ENG - US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

