



# SurgiSafe® Specialty Absorbent Floor Mats

## Intended Use


Used to help keep floors in the healthcare setting clean, dry, and slip-free particularly during fluid intense procedures.

## Key Features

- Available with and without fluid barrier backing
- Helps enable faster operating room turnaround times by quickly absorbing and retaining fluid on the floor
- Designed to help keep healthcare staff safe by reducing the chance of slips and falls due to wet floors
- Aids facilities in compliance with OSHA Bloodborne Pathogens and General Environmental Controls Standards<sup>①</sup>
- Designed to include unique features beyond those offered by the Standard Absorbent Floor Mat Line
  - Low profile
  - Offered in large sizes
- Blue mat features fluid barrier backing to help keep floor beneath dry and clean
- Yellow mat absorbs fluid from the top and bottom
- Non-sterile, intended for single-use, and not made with natural rubber latex

## SurgiSafe® Specialty Absorbent Floor Mats

Used to help keep floors in the healthcare setting clean, dry, and slip-free particularly during fluid intense procedures.



### Aspen Transfer Board

#### Technical Specifications

Part No.	Color	Absorbency	Fluid Barrier Level	Size Backing	Qty/ Case
82350	Blue	Medium	Yes	23" x 40"	50
84625	Blue	Medium	Yes	46" x 40"	25
84610	Blue	Medium	Yes	46" x 40"	10
84672	Blue	Medium	Yes	40" x 72"	10
81860-2	Yellow	Medium	No	18" x 40"	60
83630-2	Yellow	Medium	No	36" x 40"	30
83610-2	Yellow	Medium	No	36" x 40"	10
83615-2L	Yellow	Medium	No	36" x 72"	15
83625-2L	Yellow	Medium	No	36" x 72"	12
83630-M	Yellow	Medium	No	36" x 40"	15
	White	Medium	No	28" x 40"	15

SurgiSafe® is a registered trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©OSHA. Occupational Safety and Health Standards. Toxic and Hazardous Substances. Bloodborne Pathogens. 1910.1030.

©2017 Aspen Surgical, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1020114 rev 1 14-JUN-2017 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

**Phone 616-698-7100**  
**Toll-free 888-364-7004**

**[www.aspensurgical.com](http://www.aspensurgical.com)**

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**