

## WriteSite® Plus Surgical Skin Marker

### Intended Use

Intended to provide an accurate and reliable way to mark the incision site prior to surgery with Gentian Violet Ink formulated specifically for use with Chloraprep® Products.


### Key Features

- Ink is three times darker than standard Gentian Violet ink after prep with Chloraprep® products
- Ink is also compatible with Duraprep® and Betadine® prep solutions
- Ink formulated to help reduce blurring, leaving the mark clear and visible for prep
- Allows for visibility of surgical markings to support correct site surgery and the Joint Commission's Universal Protocol
- Offered in a tip style designed specifically for precision marking
- Available in regular and small barrel sizes
- Ink requires only one minute of dry time before prep for optimal results
- Virtually eliminates the need for remarking after prep
- Ink is non-tattooing and non-irritating
- Intended for single-use and not made with natural rubber latex

# WriteSite® Plus

## Surgical Skin Marker

Formulated for use with ChloroPrep® Products.



### WriteSite® Plus Standard Size Surgical Skin Marker Technical Specifications

Part No.	Description	Sterility	Qty
2700	Marker Only	Sterile	50/Box
2701	Marker, Ruler & Label Set	Sterile	50/Box
2700BN	Marker Only	Non-Sterile	100/Box
2705BN	Marker Only	Non-Sterile	50/Jar

### WriteSite® Plus Jr Size Surgical Skin Marker Technical Specifications

Part No.	Description	Sterility	Qty
2715	Marker Only	Sterile	50/Box
2716	Marker, Ruler & Label Set	Sterile	50/Box
2715BN	Marker Only	Non-Sterile	100/Jar

WriteSite® is a registered trademark of Aspen Surgical Products, Inc. ChloroPrep® is a registered trademark of CareFusion 2200, Inc. Duraprep® is a registered trademark of 3M Company. Betadine® is a registered trademark of Purdue Products L.P. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1000896 rev 11 30-Nov-2019 ENG - US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
aspensurgical.com

