



## AquaNot® Floor Suction Device

### Intended Use


Intended to help control fluids generated during wet procedures by maneuvering the AquaNot® Floor Suction Device with a push of your foot to the areas affected by fluids.

### Key Features

- The AquaNot® Device has similar characteristics as the WaterBug Quiet Floor Suction® Device, but has a lower profile
- Efficiently suction fluid waste to help reduce risks and maintain a dry operating room floor
- Suction unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- The raised design allows for easy repositioning, controlling fluids over a large floor area
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-Sterile, intended for single-use, and not made with natural rubber latex

# AquaNot® Floor Suction Device

Intended to help control fluids generated during wet procedures by maneuvering the AquaNot® Floor Suction Device with a push of your foot to the areas affected by fluids.



## AquaNot® Floor Suction Device Technical Specifications

Part No.	Description	Tubing	Qty
97020	AquaNot® Floor Suction Device	Yes	20/Case
97060	AquaNot® Floor Suction Device	Yes	60/Case
97050-NT	AquaNot® Floor Suction Device	No	50/Case

AquaNot® is a registered trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1014393 rev 3 9-Dec-2020 ENG - US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

