



## Bard-Parker® Blade Remover

### Intended Use

The Bard-Parker® Scalpel Blade Remover is intended to provide a method to remove blades from conventional scalpels.

### Key Features

- Helps safely remove up to 5 blades from scalpels without the need to touch the blade
- Helps to protect users from sharps injuries with the ability to safely unload contaminated blades into a secure channel
- Clear lid allows easy counting of the blades that have been removed
- Locking tabs prevent blades from being removed from device and clearly indicate when device is full
- Adhesive backing secures device to any surface, including drapes
- Features an audible and tactile response to ensure that blade is safely secure in device before removing scalpel handle
- Helps to comply with OSHA Bloodborne Standard 1910.1030(d)(2)(vii)(B) and AORN Recommended Practices
- Intended for single-use and not made with natural rubber latex

# Bard-Parker®

## Blade Remover

The Bard-Parker® Scalpel Blade Remover is intended to provide a method to remove blades from conventional scalpels.

### Bard-Parker® Blade Remover

#### Technical Specifications

Part No.	Description	Sterility	Color	Qty
ST-115	Blade Remover	Sterile	Red	24/Box
ST-115NS	Blade Remover	Non-sterile	Red	100/Box



Bard-Parker® is a registered trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1009574 rev 4 30-Nov-2020 ENG - US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

