



Bard-Parker® SafeSwitch® Disposable Scalpel Handle & Cover

Intended Use

Intended to be used with surgical blades for tissue separation and other procedures that require a sharp surgical blade to puncture or cut.

Key Features

- Helps provide protection from injuries that may occur before use, during steps of a procedure, or after use
- Designed for both left and right-handed surgeons
- Features a retractable shield that boasts easy one-handed activation
- Not made with natural rubber latex
- Disposable

Scalpel Handle

- Compatible with the following Bard-Parker® blade sizes: 10, 11, 15, and 15C (blades not included)
- Designed with a similar weight to existing Bard-Parker® conventional and protected blade system products
- Compatible with both Bard-Parker® stainless steel and carbon Rib-Back® blade styles

Scalpel Handle Cover


- Compatible with Bard-Parker® conventional size 3 handle (handle not included)

Bard-Parker® SafeSwitch®


Disposable Scalpel Handle & Cover

Intended to be used with surgical blades for tissue separation and other procedures that require a sharp surgical blade to puncture or cut.

Scalpel Handle



Scalpel Handle Cover



Blades sold separately.

Technical Specifications

Part No.	Description	Size	Sterility	Qty
ST-935D	Bard-Parker® SafeSwitch® Disposable Scalpel Handle	3	Sterile	8/Box, 32/Case
ST-830D	Bard-Parker® SafeSwitch® Disposable Scalpel Handle Cover	3	Sterile	8/Box, 32/Case

Bard-Parker®, Rib-Back®, and SafeSwitch® are registered trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1019934 rev 3 13-JUL-2020 ENG - US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

