

Oxboro® Corner Protectors

Intended Use


Intended to help protect instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays.

Key Features

- Aspen Surgical's corner protectors are designed to create a flexible barrier between the tray's sharp corners and the sterilization wrap
- They can be used with steam, gas, ozone, and low temperature plasma sterilization
- Reusable
- Not made with natural rubber latex

Oxboro® Corner Protectors

Intended to help protect instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays.



Specialty Disposable Light Handle Technical Specifications

Part No.	Description	Sterility	Qty
CP-80NS	Corner Protectors	Non-Sterile	80/Bag

Oxboro® is a registered trademark of Aspen Surgical. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1000885 rev 3 8-Dec-2020 ENG - US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

