



Sterion® Skin Closures

Intended Use

Intended for skin closure and closure support for sutures and staples.

Key Features

- Securely holds the catheters and other devices in place while gentle on the skin
- Available in filament reinforced and flexible styles
- No adhesive enhancers required
- Less painful to apply than suturing or stapling, providing more comfort to the patient
- Sterile, intended for single-use and not made with natural rubber latex

Filament Reinforced Skin Closures


- Reinforcing filaments offer added strength for large skin surface closures
- Filament Reinforced Skin Closure style is accepted as the standard by most healthcare professionals
- Breathable material designed to allow moisture transmission, to help reduce the chance of skin maceration
- Strips are curl-resistant

Flexible Skin Closures

- Non-woven flexible backing helps reduce the chance of skin shear that can occur when edema or flexing puts excessive pressure on the skin closure ends

Sterion® Skin Closures

Intended for skin closure and closure support for sutures and staples.



Sterion® Skin Closures – Filament Technical Specifications

Part No.	Description	mm	inches	Strips/Pouch	Pouches
751040pbx	Filament Reinforced Skin Closures	3 x 75	1/8 x 3	5	50/Box
751041pbx	Filament Reinforced Skin Closures	6 x 75	1/4 x 3	3	50/Box
751042pbx	Filament Reinforced Skin Closures	6 x 38	1/4 x 1 1/2	6	50/Box
751046pbx	Filament Reinforced Skin Closures	6 x 100	1/4 x 4	10	50/Box
751047pbx	Filament Reinforced Skin Closures	12 x 100	1/2 x 4	6	50/Box
751048pbx	Filament Reinforced Skin Closures	25.4 x 127	1 x 5	4	25/Box

Sterion® Skin Closures – Flexible Technical Specifications

Part No.	Description	mm	inches	Strips/Pouch	Pouches
753040pbx	Flexible Skin Closures	3 x 75	1/8 x 3	5	50/Box
753041pbx	Flexible Skin Closures	6 x 75	1/4 x 3	3	50/Box
753042pbx	Flexible Skin Closures	6 x 38	1/4 x 1 1/2	6	50/Box
753046pbx	Flexible Skin Closures	6 x 100	1/4 x 4	10	50/Box
753047pbx	Flexible Skin Closures	12 x 100	1/2 x 4	6	50/Box
753048pbx	Flexible Skin Closures	25 x 125	1 x 5	4	25/Box

Sterion® is a registered trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1027754 rev 2 9-Dec-2020 ENG – US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

