



## Aspen Surgical® Transparent Dressings

### Intended Use

Intended for woundcare applications, securing objects to the skin such as catheters, and as a secondary dressing.

### Key Features

- Securely holds the catheters and other devices in place while gentle on the skin
- Used to protect wound site from bacterial and viral infection
- Breathable to ensure proper moisture content in the wound bed
- Clear dressing allows for wound visualization
- Offered in two popular sizes
- Flexible and forms to body contours to provide less stress to the patient's skin during movement
- Can be worn longer than tape and gauze dressings leading to less dressing changes
- Hypoallergenic
- Available sterile or non-sterile for kit applications
- Sterile, intended for single-use, and not made with natural rubber latex

# Aspen Surgical® Transparent Dressings



For woundcare applications, securing objects to the skin such as catheters, and as a secondary dressing.

## Aspen Surgical® Transparent Dressing Technical Specifications

Part No.	Description	Sterility	Size	Qty
TDS-24	Transparent Dressing	Sterile	2.375" x 2.75"	100/Box, 4 Boxes/Case
TDS-26	Transparent Dressing	Sterile	4" x 4.5"	50/Box, 4 Boxes/Case

Aspen Surgical is a registered trademark of Aspen Surgical Products, Inc.  
Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1014595 rev 3 11-Dec-2020 ENG - US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

