



## ViraSorb® Solidifier

### Intended Use

Used to manage infectious waste fluid generated during wet procedures in surgical settings.

### Key Features

- Evenly gels canister contents in seconds without stirring or agitating
- Non-toxic, glutaraldehyde-free, and can be safely incinerated
- It does not revert back to a fluid state, helping to eliminate spillage and leakage, qualifying it as a Solid Regulated Medical Waste
- Pour nozzles are included with multi-use-bottles
- Intended for single-use and not made with natural rubber latex

# ViraSorb® Solidifier

Used to manage infectious waste fluid generated during wet procedures in surgical settings.



## ViraSorb® Solidifier Technical Specifications

Part No.	Description	Solidifies	Qty
1000	ViraSorb® Solidifier – Single-use Bottle	1,200cc	96/Case
1001	ViraSorb® Solidifier – Single-use Bottle	1,500cc	96/Case
1002	ViraSorb® Solidifier – Single-use Bottle	2,000cc	36/Case
1003	ViraSorb® Solidifier – Single-use Bottle	3,000cc	36/Case
1005	ViraSorb® Solidifier – Single-use Pouch	500cc	200/Case
1010	ViraSorb® Solidifier – Multi-use Bottle	10,000cc	12/Case
1020	ViraSorb® Solidifier – Multi-use Bottle	20,000cc	12/Case

ViraSorb® is a registered trademark of Aspen Surgical. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1014679 rev 3 9-Dec-2020 ENG – US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

