

Personal Protective Equipment

Accessories

 **aspen** surgical®



Aspen Surgical strives for unmatched quality and unparalleled dependability with brands like Bard-Parker®, Colby™, Precept®, Protek®, Stork®, SurgiSafe®, WriteSite®, and more. Aspen Surgical is a leading provider of surgical essentials, PPE, and ultrasound imaging accessories. We got here by providing solutions that are as effective as they are trusted. With decades of expertise in the operating rooms, we are not always recognized by name, but our reputation is always at hand. And although we have been supporting medical teams for many years, we never stop evolving, looking ahead, and expanding our reach.

Learn more about what Aspen Surgical has to offer in the pages of this catalog, or online at www.aspensurgical.com



Table of Contents

Headwear... 6-8

- Bouffant, Skull Cap and Surgeon Caps pg... 6
- Surgical Hoods and Technique Cap pg... 7
- Tyvek® Hood pg... 8
- Beard Cover pg... 8

Footwear... 9

- Wave™ Shoe Covers pg... 9
- Tyvek® Shoe and Boot Covers pg...9

Lab Coats and Jackets... 10

- Lab Jacket pg... 10
- Lab Coat pg... 10

Coveralls... 11

- SMS Coverall pg... 11
- Tyvek® Coverall pg... 11

PPE Apparel... 12

- Patient Gown pg... 12
- Apron pg... 12

Miscellaneous Accessories... 13-14

- Tyvek® Sleeve pg... 13
- Infant Accessories pg... 13
- Spill Pad pg... 14

Cold Packs... 14

- Premium Cold Packs pg... 14
- Cold Pack pg... 14



Aspen Surgical Products, Inc. Headquarters, Caledonia, Michigan

DuPont™ and Tyvek® are trademarks or registered trademarks of DuPont or its affiliates.

The following brands are trademarks or registered trademarks of Aspen Surgical Products, Inc.: Precept® and Wave™

© 2021 Aspen Surgical Products, Inc. All rights reserved.

Aspen Surgical’s line of headwear is designed to meet the various needs of healthcare personnel. Headwear is available in a variety of styles, including bouffants, surgical caps, surgeon hoods, and beard covers.

Headwear



- Bouffant Cap**
Product Details
- Designed for user comfort
 - Elastic headband
 - Made with polypropylene material
 - Made in North America
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
22133	X-Large	White	24" Bouffant Cap with Elastic Headband	Polypropylene	125/Box	500/Case



- Skull Cap**
Product Details
- Designed for user comfort
 - Soft absorbent headband
 - Elastic back keeps hair in place while providing a flexible fit
 - Made with lightweight SMS material
 - Made in North America
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
1014P	Universal	Blue	Skull Cap with Soft Absorbent Headband	SMS	125/Box	500/Case



- Surgeon Cap**
Product Details
- Made with absorbent Kaycel fabric
 - Designed for user comfort
 - Made in North America
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
5200	Universal	Blue	Surgeon Cap	Kaycel	125/Box	500/Case



- Surgical Hood**
Product Details
- Soft SMS material
 - Tie closure
 - Increased facial coverage
 - Made in North America
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
1011P	Universal	Blue	Surgical Hood with Tie Closure	SMS	100/Box	200/Case



- Legionnaire Surgical Hood**
Product Details
- Soft SMS material
 - Tape closure
 - Increased facial coverage
 - Made in North America
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
5246	Universal	Blue	Legionnaire Surgical Hood with Tape Closure	SMS	100/Box	200/Case



- Technique Cap**
Product Details
- Soft SMS material
 - Tape closure
 - Universal sized
 - Made in North America
 - 100/Box, 300/Case
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
2210	Universal	Blue	Technique Cap with Tape Closure	SMS	100/Box	300/Case

Headwear - Tyvek® Hood (Clean Room)



- Tyvek® Hood**
Product Details
- Full face with snaps
 - Large
 - Made in North America
 - 100/Case Bulk
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
35-60230	Large	White	Tyvek® Hood with Snap Closure	Tyvek®	Bulk	100/Case

Headwear - Beard Cover



- Beard Cover**
Product Details
- Lightweight polypropylene material
 - Soft comfortable elastic headband
 - Excellent breathability
 - Made in North America
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
35-3015	Universal	White	Beard Cover with Elastic Headband	Polypropylene	500/Box	500/Case

Aspen offers shoe and boot covers with anti-skid embossed polyethylene sole, constructed of material such as polypropylene or Tyvek® by DuPont, and the convenience of dispenser box or bulk packaging.

Shoe Covers



- Wave™ Shoe Covers**
- Anti-skid embossed polyethylene sole
 - Triple-elastic runs that allow a smooth wrap around the foot to prevent tripping
 - Available in Easy-Open bags, Easy Pull-Out or Econo dispensing boxes
 - Made in North America
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
SHS-0202-MF	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	3 Easy-Open Bags	300/Case
SHS-0202-CB	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	6 Easy-Pull-Out Dispenser Boxes	600/Case
SHS-6202-MF	X-Large	Blue	Shoe cover anti-skid sole	Polypropylene	3 Easy-Open Bags	300/Case
SHS-6202-CB	X-Large	Blue	Shoe cover anti-skid sole	Polypropylene	2 Econo Dispenser Boxes	400/Case

Tyvek® Shoe and Boot Covers




- Tyvek® Shoe and Boot Covers**
Product Details
- Anti-skid embossed polyethylene sole
 - Large, X-Large
 - Made in North America
 - 100/Case Bulk
 - Not made with natural rubber latex

Part No.	Size	Color	Style	Material	Packaging	Qty
35-60420L	Large	White	Shoe cover anti-skid sole	Tyvek®	Bulk	100/Case
35-60420XL	X-Large	White	Shoe cover anti-skid sole	Tyvek®	Bulk	100/Case
35-60510L	Large	White	Boot cover 18" tall, anti-skid sole	Tyvek®	Bulk	100/Case
35-60510XL	X-Large	White	Boot cover 18" tall, anti-skid sole	Tyvek®	Bulk	100/Case

Lab Coat and Jacket
Product Details

- Soft knit cuffs and collar provide maximum comfort while enhancing fit
- Snap closures and 3 pockets provide efficiency and convenience
- Special finish provides durable tear resistant material that allows for ease in movement
- Sizing from petite to double extra large
- Intended to be used in environments where low fluid exposure is a risk
- Made in North America
- Not made with natural rubber latex

Sizing Chart — Approximation Only




Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
1239P	Petite	39½"	48"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239M	Medium	40½"	50"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239L	Large	41½"	52"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239XL	X-Large	42½"	54"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239XXL	XX-Large	43½"	56"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725P	Petite	29"	48"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725M	Medium	29¾"	50"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725L	Large	31¼"	52"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725XL	X-Large	33"	54"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725XXL	XX-Large	31¾"	58"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case

Coveralls, SMS and Tyvek®
Product Details

- Coveralls are available in a variety of styles and material such as multi-layer SMS or Tyvek®
- Special finish SMS material provides durable tear resistance that allows for ease in movement
- Low-linting, anti-static white Tyvek® material is noted for its strength and resistance to tears and punctures,
- Generous sizing
- Full collar with elastic or open wrist and ankles
- Made in North America
- Not made with natural rubber latex

Sizing Chart — Approximation Only



Part No.	Size	Color	Style	Material	Packaging	Qty
1547M	M-XXXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547L	M-XXXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547XL	M-XXXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547XXL	M-XXXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547XXXL	M-XXXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1548M	M-XXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1548L	M-XXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1548XL	M-XXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1548XXL	M-XXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
35-0429XXL	XX-Large	White	Full collar with elastic wrists and ankles	Tyvek®	25/Bag	25/Case
35-1412XL	X-Large	White	Full collar with elastic wrists and ankles	Tyvek®	25/Bag	25/Case



Patient Gown, Full Coverage

Product Details

- To maintain patient dignity the gown has 3/4 sleeves with ties at neck and waist
- Made with multi-layer SMS material
- Good opacity and tear resistance
- Disposable patient gowns help provide infection prevention and help eliminate the concern of missing buttons or ties
- Made in North America
- Not made with natural rubber latex

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
0456P	Universal	47"	57"	Dark Blue	Secure Neck and Waist Ties	SMS	100/Bag	100/Case

Apron



Apron, Blue Film Over-the-Head Open Back

Product Details

- Heavyweight blue film material
- Over-the-head donning
- Secure waist ties
- Die cut
- Made in North America
- Not made with natural rubber latex

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
8571	Universal	59"	43"	Blue	OTH Open Back	Blue Film	20/Bag	20/Case



Tyvek® Protective Sleeve

Product Details

- Lightweight Tyvek® material provides excellent protection against fluids and bacteria
- Tyvek® is noted for its strength and resistance to tears and punctures
- Low-linting, anti-static
- Made in North America
- Not made with natural rubber latex

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
35-2502	Universal	18"	NA	White	Sleeve with elastic	Lightweight Tyvek®	200/Bag	200/Case

Newborn Accessories




Newborn Accessories

Product Details

- Disposable newborn accessories are designed with the comfort of infants in mind
- Materials are soft, yet durable, to help provide a nurturing environment
- Made in North America
- Not made with natural rubber latex

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
5180	Newborn	NA	NA	Pink and Blue Print	Infant Shirt with Tape Closure	Spunlace	100/Box	400/Case
5431	NA	19"	15"	Blue	Infant Scale Line	SMS	250/Bag	1000/Case
60112	NA	30"	36"	White	Receiving Blanket	Rayon/ Polypropylene Blend	100/Bag	100/Case



Spill Pad

Product Details

- Intended to be used for the recovery of spilled liquids and clean-up
- Absorbent material
- Fluid impervious on one side
- Made in North America
- Not made with natural rubber latex

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
99226	NA	13"	16"	Blue	Spill Pad	Absorbent blend spunlace laminated to polypropylene blend	200/Bag	200/Case

Cold Packs

Cold Packs

Product Details

- Available in three sizes, premium cold packs, large and small, and universal sized cold pack
- Wide mouth with opening aids that allow for extremely easy filling
- Clip sleeve provides a secure holding area for the clip
- Secure clip closure
- Sturdy and strong ties allow for secure attachment
- Three-layer design, soft outer layer for comfort, absorbent middle layer and an impervious inner plastic bladder which helps provide even cold therapy
- Not made with natural rubber latex

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
16100	Large	13"	8"	White	Premium Cold Pack	Rayon/Polyester	25/Box	100/Case
16101	Small	13"	6"	White	Premium Cold Pack	Rayon/Polyester	25/Box	100/Case
16101	Universal	13"	7½"	White	Cold Pack	Rayon/Polyester	25/Box	100/Case

Notes:

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

DuPont™ and Tyvek® are trademarks or registered trademarks of DuPont or its affiliates.

Precept® and Wave™ are trademarks or registered trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2022 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1057537 rev 2 27-JUN-2022 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

USA 616.698.7100
Canada 888.364.7004
www.aspensurgical.com

