



## 65-3323 Tape Fog-Shield® Surgical Mask

### Intended Use

Intended to be used in environments where high levels of fluid exposure is a risk.



### Complies with ASTM F2100

Fluid Resistance (mmHg): 160  
BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 98\%$   
PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 98\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 5$   
Flame Spread: Class 1

### Key Features

- Complies with ASTM F2100 Level 3 requirements
- Fog-Shield® anti-fog tape inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Strong polyester ties provide a secure and comfortable fit
- Made in North America
- Not made with natural rubber latex

# Tape Fog-Shield® Surgical Mask

Intended to be used in environments where high levels of fluid exposure is a risk.



## Technical Specifications

Part No.	65-3323
Description	Tape Fog-Shield® Surgical Mask
Color	Green Diamond
Fluid Resistance	160 mm Hg
Filtration – BFE	≥99.8
Filtration – PFE	≥99.0
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Vertical Polyester Ties
Mask Construction	Sewn
Qty	50/Box, 5/Boxes, 250/Case

Precept® and Fog-Shield® are registered trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and model.

©2023 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1110421 rev 1 28-MAR-2023 ENG – US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

