

# Orthopedic Surgery Disposables

High Quality, Low Cost,  
American Made

 **aspen** surgical®



We're Aspen Surgical, a supplier of surgical essentials. From the precision of the finest scalpel with our Bard-Parker® brand to our Precept® PPE and Protek™ ultrasound probe covers that help fight preventable infections and cross contamination. And our WriteSite® Plus skin markers that provide an accurate and reliable way to mark the incision site before surgery.

From one trusted supplier, you get a portfolio of essential, physician preferred medical disposables, sold through an integrated channel. With strong distributor relationships and a dedicated field service team, we ensure we stay grounded in our customer's needs every step of the way.

We help provide an easier set-up, clean up, and quicker turnover in the O.R. for surgical teams and facilities so you can focus on what matters the most – your patients.

From facilities across North America, our dedicated team designs, manufactures, and delivers high-quality products to healthcare facilities around the world.

Every day we strive to make a difference. Every day we hold ourselves to a higher standard, so whenever and wherever you need us, we'll be there.

Learn more about all that Aspen Surgical has to offer in the pages of this catalog, or online at [www.aspensurgical.com](http://www.aspensurgical.com)







Aspen Surgical Products, Inc. Corporate Headquarters, Caledonia, Michigan

The following brands are registered trademarks of Aspen Surgical Products, Inc:  
CinchMat®, ErgoSupport®, EvacMat®, Hoffee®, Precept®, Strako®, SurgiSafe®, WaterBug Quiet Floor Suction®, and WriteSite®.

The following brands are registered trademarks of Aspen Surgical Products, Inc by Colby Manufacturing:  
AquaNot®, ViraSorb®, and WaterBoom®.

The following brands are trademarks of Aspen Surgical Products, Inc:  
Colby™, Hermanator™, Protek™

Trumpf Medical™ is a trademark of Trumpf Medizin Systeme GmbH & Co. KG and iLED® is a registered trademark of Trumpf Medizin Systeme GmbH & Co. KG.  
TruLight™ is a trademark of Trumpf GmbH + Co. KG.

Please note, not all products are available in all regions.

©2023 Aspen Surgical Products, Inc. All rights reserved.

# Table of Contents

Minimally Invasive... 7-21

Hoffee® Blue Lasso Suture Retriever . . . . .7

Green Genie Suture Passer. . . . .7

Orthopedic Patient Positioning... 8-21

Erin Beach Chair Face Mask . . . . .8

Erin Beach Chair Face Mask Replacement Straps . . . . .8

Universal Head Restraint . . . . .8

Indy Beach Chair Face Mask. . . . .9

Bethesda Beach Chair Face Mask . . . . .9

Beach Chair Head and Chin Straps . . . . .9

Marco Shoulder Stabilization Kit. . . . .10

Berkette Beach Chair Limb Positioner Kit . . . . .10

Jessie Arm Suspension Kit . . . . .11

Norton Arm Restraint Kit . . . . .11

Dr Don’s® Arm Canoe . . . . .12

Dr Don’s® Arm Canoe-V . . . . .12

Dr Don’s® Arm Suspension Kit . . . . .13

Dr Don’s® Arm Suspension Kit – Montana . . . . .13

Dr Don’s® Arm Suspension Kit – Northern . . . . .13

Carson 8" Perineal Post Pad . . . . .14

Patient Care Boot Liner Pad. . . . .14

Perineal Post Cover . . . . .14

Lewiston Positioner Kit . . . . .15

Lindy Traction Boot Liner . . . . .15

Lateral Positioner Pad . . . . .15

Racine Pelvic Positioning Pad . . . . .16

Hermanator® Knee Positioning Pad. . . . .16

XL Hermanator® Pink Knee Positioning Pad . . . . .17

Strako® Knee Positioning Pad . . . . .17

TLC Distraction Pad Set . . . . .18

Leg Comfort Pad . . . . .18

Smoot Ankle Distraction Strap . . . . .19

Thigh Holder Pad . . . . .19

Parent Post Knee Pad . . . . .19

Prone Face Pillow . . . . .20

Sedona Arm Positioner Set . . . . .20

MD Traction Rope . . . . .21

Brachytherapy Mattress . . . . .21

Fluid Management & Floor Safety... 22-30

SurgiSafe® Specialty Absorbent Floor Mats . . . . .23

SurgiSafe® Standard Absorbent Floor Mats . . . . .24

SurgiSafe® EvacMat® Absorbent Floor Mat. . . . .25

SurgiSafe® CinchMat® Absorbent Floor Mat. . . . .25

ErgoSupport® Anti-Fatigue Mat . . . . .26

Entrance Mat (TacTec) . . . . .26

ErgoSupport® Anti-Fatigue Step Stool Mat . . . . .27

ErgoSupport® Anti-Fatigue Step Stool . . . . .27

WaterBug Quiet Floor Suction® Device . . . . .28

WaterBoom® Floor Suction Strip . . . . .28

AquaNot® Floor Suction Device. . . . .29

ViraSorb® Solidifier . . . . .29

AspenBasin® Fluid Waste Container . . . . .30

Products by Original Manufacturer... 34-35



Suture Retrievers

The North American-made suture retrievers below provide high quality for the clinical staff, while at the same time offering cost savings to the medical facility. The knurled aluminum handles and distinctive coloring for ease of use and comfort in-hand. These suture retrievers are designed for single patient use. Intended to aid in the passing of sutures through bone and soft tissue during orthopedic surgery.



- Hoffee® Blue Lasso Suture Retriever**  
**Product Details**
- Comparable to Smith & Nephew® Hewson® Suture Retriever 7111-1579 and Arthrex® Transosseous Suture Passer AR-1000
  - Hewson Suture Retriever replacement
  - Aluminum knurled blue handle is weighted for ease of use
  - Blue suture loop for increased visibility
  - 10" length is suitable for most laparoscopic procedures

\*Please inquire about Non-sterile suture retrievers for custom packs

Part No.	Sterility	Qty
710000	Sterile, intended for single-use	5 Eaches/Box



- Green Genie Suture Passer**  
**Product Details**
- Comparable to SWANSON™ Suture Passer
  - Aluminum knurled green handle is weighted for ease of use
  - Wire suture loop
  - 6.5" length
  - Made in North America

Part No.	Sterility	Qty
710001	Sterile, intended for single-use	5 Eaches/Box



Beach Chair Accessories

The North American-made head restraints are used with the following beach chairs: Allen®, Arthrex®, SchureMed™, Skytron®, Smith & Nephew® and Steris®. Each surgical positioning foam product is designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



Erin Beach Chair Face Mask

Product Details

- For use with Smith & Nephew® T-MAX™ Beach Chair 7210996
- Comparable to Smith & Nephew® T-MAX™ Disposable Face Mask 7210559
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711001	Non-sterile, intended for single-use	10/Box



Erin Beach Chair Face Mask Replacement Straps

Product Details

- For use with Smith & Nephew® T-MAX™ Beach Chair 7210996
- Comparable to Smith & Nephew® T-MAX™ Replacement Strap 7210558
- Compatible with Erin Beach Chair Face Mask (711001)
- Reusable

Part No.	Sterility	Qty
711002	Non-sterile	1 Pair/Box



Universal Head Restraint

Product Details

- For use with Allen® Beach Chair Positioner and Arthrex® Beach Chair Positioner
- Comparable to Allen® Universal™ Head Restraint A-90023 and Arthrex® Universal Head Restraint AR-1627-06
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711004	Non-sterile, intended for single-use	10/Box



Indy Beach Chair Face Mask

Product Details

- For use with Allen® Flexiform™ Head Rest A-90017
- Comparable to Allen® Flexiform™ Head Rest Disposable A-90019
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711003	Non-sterile, intended for single-use	10/Box



Bethesda Beach Chair Face Mask

Product Details

- For use with SchureMed™ Beach Chair 800-0004, Stryker® Shoulder Chair OZ 7519475 and OZ 7523075, and Skytron® Beach Chair
- Comparable to SchureMed™ Disposable Face Mask 800-0164, Stryker® Disposable Mask OZ 7547575, and Skytron® Disposable Face Mask 3-040-33
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711000	Non-sterile, intended for single-use	10/Box



Beach Chair Head and Chin Straps

Product Details

- For use with Skytron® Beach Chair, SchureMed™ Beach Chair, Steris® Beach Chair, Stryker® Shoulder Chair OZ 7519475 and OZ 7523075
- Comparable to Skytron® Beach Chair Head & Chin Straps 3-040-22, SchureMed™ Disposable Head & Chin Strap 800-0039, Steris® Head Restraint Straps BF039, and Stryker® Chin Strap OZ 7518375
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711005	Non-sterile, intended for single-use	20/Box

Shoulder Positioner Kits

The North American-made sterile kits provide high quality for the clinical staff, while at the same time offering cost savings to the medical facility. For example, the lateral kits feature a snap clip to more securely attach the product to the frame and cost less than the brand name. Each kit is designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.




Marco Shoulder Stabilization Kit

Product Details

- For use with Smith & Nephew® SPIDER™ Limb Positioner 7210569
- Comparable to Smith & Nephew® SPIDER™ Shoulder Stabilization Kit 7210573
- Includes Drape, Arm & Hand Cuffs, Hand Grip, and 4" x 5yd 3M® Coban™
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711100	Sterile, intended for single-use	12/Box



Berkette Beach Chair Limb Positioner Kit

Product Details

- For use with Arthrex® Trimano® Shoulder Positioner AR-1640
- Comparable to Arthrex® Trimano® Beach Chair Shoulder Kit AR-1644
- Includes Foam Sleeve, Hand Grip, Drape, and 4" x 5yd 3M® Coban™
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711200	Sterile, intended for single-use	6/Box



Jessie Arm Suspension Kit

Product Details

- For Use With McConnell Limb Positioner
- Comparable To McConnell Arm Support 12-401
- Includes Arm Wrap, Placement Strap, and Hand Hexapod
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711300	Sterile, intended for single-use	10/Box



Norton Arm Restraint Kit

Product Details

- For Use With Allen® Beach Chair A-92000
- Comparable To Allen® Arm Positioner Restraint A-92001
- Includes Foam Arm Drape
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711400	Sterile, intended for single-use	10/Box



Dr Don's® Arm Canoe

Product Details

- For Use With Arthrex® Lateral Decubitus Shoulder Traction Tower AR-1630 and Arthrex® 3-Point Shoulder Distraction System AR-1600M
- Comparable To Arthrex® STAR™ Sleeve AR-1606
- Includes Convuluted Foam Sleeve, Arm Support Wrap, and 4" x 5yd 3M® Coban™
- Soft, non-abrasive foam for patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711500	Sterile, intended for single-use	5/Box



Dr Don's® Arm Canoe-V

Product Details

- For Use With Arthrex® Lateral Decubitus Shoulder Traction Tower AR-1630 and Arthrex® 3-Point Shoulder Distraction System AR-1600M
- Comparable To Arthrex® STAR™ Sleeve AR-1606V
- Includes Convuluted Foam Sleeve and Arm Support Wrap
- Soft, non-abrasive foam for patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711501	Sterile, intended for single-use	5/Box



Dr Don's® Arm Suspension Kit

Product Details

- Comparable To Instramed® Upper Extremity Drape IN-9300CLF and DeRoyal® Shoulder Suspension Kit 32-525
- Includes fluid resistant sleeve, rope and snap hook, and 4" x 5yd 3M® Coban™
- Fluid resistant stockinette allows for secure patient positioning
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711502	Sterile, intended for single-use	10/Box



Dr Don's® Arm Suspension Kit – Montana

Product Details

- Comparable To Smith & Nephew® Shoulder Suspension Kit 72200195 and Zimmer Biomet® Shoulder Suspension Kit 905026
- Includes Reinforced sleeve, rope and snap hook and 4" x 5yd 3M® Coban™
- Reinforced, soft cloth lined stockinette allows for secure patient positioning and patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711503	Sterile, intended for single-use	10/Box



Dr Don's® Arm Suspension Kit – Northern

Product Details

- Comparable To Arthrex® Lateral Arm Traction Sleeve AR-1635
- Includes Jersey knit traction sleeve, grip lining, rope and snap hook, 4" x 5yd 3M® Coban™ and foam hand grip
- Jersey knit stockinette with interior grip lining allows for secure patient positioning
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711504	Sterile, intended for single-use	10/Box



Hip Positioners

The North American-made positioning foam products below are compatible with Innovative Medical Products® (IMP®), Mizuho OSI® and Smith & Nephew® hip positioning devices. These positioning foam products are designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



Carson 8" Perineal Post Pad

Product Details

- For Use With Smith & Nephew® Supine Hip Positioning System 72200624 and Mizuho OSI® Hana®, ProFx®, and Hana SSXT® Tables
- Comparable To Smith & Nephew® Supine Perineal Post Pad 72200634
- Soft, non-abrasive foam for patient comfort
- Outside diameter: 9", Height: 11" Core Diameter: 1.5"

Part No.	Sterility	Qty
712200	Non-sterile, intended for single-use	1/Box



Patient Care Boot Liner Pad

Product Details

- For Use With Smith & Nephew® Supine Hip Positioning System 72200624 and Smith & Nephew® SPIDER™ Boot 72201425
- Comparable To Smith & Nephew® Boot Insert Pads 72200635 and Smith & Nephew® SPIDER® Boot Liner 72201426
- Soft, non-abrasive foam for patient comfort
- IFU printed on product for quick reference

Part No.	Sterility	Qty
712300	Non-sterile, intended for single-use	1 Pair/Box



Perineal Post Cover

Product Details

- For Use With Mizuho OSI® Hana®, ProFx®, and Hana SSXT® Tables
- Comparable To Mizuho OSI® Perineal Post Cover 6855-13
- Soft, non-abrasive foam for patient comfort
- Dimensions: 11.5" x 5.5"

Part No.	Sterility	Qty
712400	Non-sterile, intended for single-use	12/Box



Lewiston Positioner Kit

Product Details

- For Use With Arthrex® Hip Distraction System AR-6529S
- Comparable To Arthrex® Disposable Hip Kit AR-6529-11
- Includes Sueded perineal post pad cover and 1 pair of traction boot pads
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
712100	Non-sterile, intended for single-use	1/Box



Lindy Traction Boot Liner

Product Details

- For Use With Allen® Traction Boot A-93005 and A-93000 and Steris® OT 1000 Traction Boot
- Comparable To Allen® Traction Boot Disposable A-93001 and STERIS® Traction Boot Liner BF419
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
712900	Non-sterile, intended for single-use	6/Box




Lateral Positioner Pad

Product Details

- For Use With IMP® Universal Lateral Positioner® 1092
- Comparable To IMP® Universal Lateral Positioner® Economy Patient Protective Pad Set 1092-P
- Soft, non-abrasive foam for patient comfort
- Includes anterior support pads

Part No.	Sterility	Qty
712600	Non-sterile, intended for single-use	15/Box





Racine Pelvic Positioning Pad

Product Details

- For Use With IMP® McGuire™ Hip Positioner 1084
- Comparable To IMP® McGuire™ Hip Pads 1084-P
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
712800	Non-sterile, intended for single-use	15/Box

Knee & Ankle Positioners

Each product below provides high quality for the clinical staff, while at the same time offering cost savings to the medical facility. These sterile kits and positioning foam products are designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



Hermanator® Knee Positioning Pad

Product Details

- For Use With IMP® De Mayo Knee Positioner® 803 AND 803-S, and Match Grade Ladderloc Knee Positioner, SchureMed™ TKR Knee Positioner, and Innomed® Robb Knee Positioner
- Comparable To IMP® De Mayo Knee Positioner Pad 294-P-10 and 803-GP-10, Zimmer® Alvarado Disposable Liner 00-1320-015-10, Innomed® Sterile Protective Pad 2629-00, and SchureMed™ TKR Disposable Liner 800-0146
- Helps decrease device related mechanical pressure injury, especially to the high risk areas of the posterior, medial, and lateral malleolus
- Soft, high-density memory foam provides enhanced patient safety
- Not made of natural rubber latex

Part No.	Sterility	Qty
713100	Sterile, intended for single-use	10/Box




XL Hermanator® Pink Knee Positioning Pad

Product Details

- For use with IMP® De Mayo Knee Positioner® 803 AND 803-S and Match Grade Ladderloc Knee Positioner
- Comparable to IMP® De Mayo Knee Positioner® 803-GP-10
- Helps decrease device related mechanical pressure injury, especially to the high risk areas of the posterior, medial, and later malleolus
- Soft, high-density memory foam provides enhanced patient safety
- Not made of natural rubber latex

Part No.	Sterility	Qty
713300	Sterile, intended for single-use	10/Box



Strako® Knee Positioning Pad

Product Details

- For use with the Stryker® Mako® Leg Positioner
- Comparable to Stryker® Mako® Leg Positioner Disposable Kit
- Secures patient's foot into the leg holder boot to help stabilize the knee during surgery
- Helps decrease device related mechanical pressure injury
- Sterile, intended for single-use and not made with natural rubber latex

Part No.	Sterility	Qty
713800	Sterile, intended for single-use	6/Box



TLC Distraction Pad Set

Product Details

- For Use With IMP® De Mayo Universal Distractor® 907
- Comparable To IMP® De Mayo Pressure Protector Pad® Set 907-GP-10
- Includes Memory pad, 4" x 5yd 3M® Coban™ and Knee Distraction Guard
- Semi circle memory pad prevents device from rolling into the popliteal area while providing a safe view into the posterior capsule
- Helps protect the posterior thigh from mechanical pressure injury
- Soft, high-density memory foam provides enhanced patient safety
- Not made of natural rubber latex

Part No.	Sterility	Qty
713200	Sterile, intended for single-use	10/Box



Leg Comfort Pad

Product Details

- For Use With Allen® Leg Holder R-734-33, Smith & Nephew® Acufex Knee Holder, SchureMed™ Schure Grip Leg Holder 800-0065, and Stryker® Arthroscopy Leg Holder OZ 7511975
- Comparable To Allen® Leg Holder Disposable R-S4133-V, Smith & Nephew® Disposable Leg Holder Pads 3757, SchureMed™ Schure Grip Leg Holder Pads 508-0103, and Stryker® Disposable Pad OZ 7512075
- Inner foam is soft and non-abrasive for patient comfort
- Outer foam is a dense material that helps provide appropriate clamping pressure to stabilize the leg

Part No.	Sterility	Qty
713600	Non-sterile, intended for single-use	12/Box



Smoot Ankle Distraction Strap

Product Details

- For Use With Smith & Nephew® GUHL™ Positioner and Arthrex® Ankle Distractor
- Comparable To Smith & Nephew® GUHL™ Ankle Distraction Strap 014407 and Arthrex® Ankle Distraction Strap AR-1712
- Soft, non-abrasive foam where ankle meets strap for patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure

Part No.	Sterility	Qty
713400	Sterile, intended for single-use	10/Box




Thigh Holder Pad

Product Details

- For Use With Smith & Nephew® FERKEL™ Thigh System 72201812
- Comparable To Smith & Nephew® FERKEL™ Thigh System 72201814
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
713500	Non-sterile, intended for single-use	6/Box



Parent Post Knee Pad

Product Details

- For Use With Arthrex® PivotPost® AR-1560
- Comparable To Arthrex® PivotPost® AR-1560-01
- Soft, non-abrasive foam for patient comfort
- Velcro® straps allow for efficient placement on PivotPost®

Part No.	Sterility	Qty
712500	Non-sterile, intended for single-use	5/Box

Spine Positioning Foam

The North American-made positioning products are designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.




Prone Face Pillow

Product Details

- Comparable To Allen® Foam Head Block A-70862
- Consists of two specifically designed foam layers: bottom layer has supportive foam with pre-cut dual intubation slots with the top layer consisting of soft, non-abrasive foam designed for patient comfort
- Can be used with the C-Flex® Head Positioning System, Allen® Foam Block Stand, and the Allen® Advance Table Head Plate or any flat surface
- 5" Foam
- Memory Foam Topping
- Compressed Packaging

Part No.	Sterility	Qty
714100	Non-sterile, intended for single-use	12/Box



Sedona Arm Positioner Set

Product Details

- Comparable To Mizuho OSI® Laminectomy Arm Cradle 5823-10
- Soft, non-abrasive foam pads the entire length of the patient's arm and has notches to allow for positioning of the patient's arm
- Dimensions: 24" x 5" x 3"
- Can be used in prone or supine position

Part No.	Sterility	Qty
714200	Non-sterile, intended for single-use	6 Pairs/Box

Orthopedic Accessory Products

The North American-made sterile kits and positioning foam products are designed for specialty patient positioning needs and single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.




\*Please inquire about Non-sterile MD Traction Rope for custom packs

MD Traction Rope

Product Details

- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure Snap Hook
- Rope Length Approx 89"

Part No.	Sterility	Qty
715000	Sterile, intended for single-use	50/Box



Brachytherapy Mattress

Product Details

- Soft, non-abrasive foam for patient comfort
- 30" x 79" x 5"
- 18" x 6" Opening

Part No.	Sterility	Qty
715100	Non-sterile, intended for single-use	1/Box






Colby™ Fluid Management Products

Our Colby™ line provides a convenient and effective means of managing infectious waste fluid generated during wet procedures. Our products help to reduce the likelihood of injuries associated with wet floors. They aid in protecting healthcare professionals and housekeeping personnel from contact with contaminated substances and help support your efforts to comply with OSHA’s Bloodborne Pathogens Standard. Our products do not change the contents of medical waste, but rather help your team to more effectively contain it.

SurgiSafe® Absorbent Floor Mats

- Help keep floors clean, dry, and slip-free, even during fluid intense procedures
- Contains fluid in one place to help make operating room turnover time faster and more efficient
- Absorb irrigates, blood, and bodily fluids through unique wicking, absorbency, and fluid retention properties, while conforming to uneven floor surfaces without bunching
- Help protect your staff from slips, trips, and falls in the operating room
- Help reduce staff exposure to infectious waste
- Helps eliminate using blankets and towels which can cause tripping hazards
- Non-sterile, intended for single-use, and not made with natural rubber latex



SurgiSafe® Specialty Absorbent Floor Mats

Product Details

- Designed for strength and durability
- Features low profile design, which allows for foot pedals or stools to roll with ease
- Offered in large sizes to accommodate surgical procedures where more than one person needs to stand on the mat
- Blue mat features fluid barrier backing to help keep floor beneath dry and clean
- Yellow mat absorbs fluid from the top and bottom

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
82350	Blue	Medium	Yes	21" x 44"	50/Case
84625	Blue	Medium	Yes	32" x 46"	25/Case
84610	Blue	Medium	Yes	32" x 46"	10/Case
84672	Blue	Medium	Yes	32" x 72"	10/Case
81860-2	Yellow	Medium	No	16" x 45"	60/Case
83630-2	Yellow	Medium	No	32" x 45"	30/Case
83610-2	Yellow	Medium	No	32" x 45"	10/Case
83615-2L	Yellow	Medium	No	32" x 80"	15/Case
83625-2L	Yellow	Medium	No	32" x 80"	12/Case
83630-M	Yellow	Medium	No	32" x 45"	15/Case
	White	Medium	No	28" x 40"	15/Case





**SurgiSafe® Standard Absorbent Floor Mats**  
 Product Details

- Made with a minimum of 60% pre-consumer recycled material
- Available in three absorbency levels as well as with and without fluid barrier backing for convenience
- Pink – High level of fluid absorbency, will absorb 16x’s its dry weight
- White – Medium level of fluid absorbency, will absorb 14x’s its dry weight
- Green – Low level of fluid absorbency, will absorb 11x’s its dry weight

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
83430	Pink	High	No	28" x 40"	30/Case
83472	Pink	High	No	28" x 72"	30/Case
83330	Pink	High	Yes	28" x 40"	30/Case
83372	Pink	High	Yes	28" x 72"	30/Case
81860	White	Medium	No	18" x 40"	60/Case
818100	White	Medium	No	18" x 21"	100/Case
81850	White	Medium	No	18" x 21"	50/Case
83610	White	Medium	No	28" x 40"	10/Case
83630	White	Medium	No	28" x 40"	30/Case
83625-L	White	Medium	No	28" x 72"	12/Case
83672	White	Medium	No	28" x 72"	30/Case
83530	White	Medium	Yes	28" x 40"	30/Case
83572	White	Medium	Yes	28" x 72"	30/Case
84030	Green	Low	No	28" x 40"	30/Case
83030	Green	Low	Yes	28" x 40"	30/Case
83072	Green	Low	Yes	28" x 72"	30/Case
830125	Green	Low	Yes	28" x 125"	24/Case



**SurgiSafe® EvacMat® Absorbent Floor Mat**  
 Product Details

- Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor
- Extremely quiet operation helps to decrease distraction of the O.R. staff- as quiet as an indoor conversation
- The top layer of the mat contains wicking to help prevent pooling
- Works well with wall suction
- Simple to use and clean up – red bag friendly
- Low profile allows for the O.R. staff to roll equipment over the mat without loss of suction
- Large enough for two people
- Provides cushion to help reduce fatigue during long procedures
- Aids in keeping the O.R. floor dry to reduce the risk of injuries elated to wet floors
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
95005	Yellow	Medium	Yes	36" x 48"	5/Case



**SurgiSafe® CinchMat® Absorbent Floor Mat**  
 Product Details

- Combines the effective absorbency of the SurgiSafe® mat with an easier way to handle the mat during disposal
- Pull cord design helps make clean-up and O.R. turnaround time simple and fast
- The top layer has wicking abilities that quickly absorbs irrigants, blood, and bodily fluids from the O.R. floor to help prevent pooling and splashing
- Skid-resistant fluid barrier backing helps ensure that the floor beneath the mat remains dry and prevents staining
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
84625-C	Blue	Medium	Yes	32" x 46"	10/Case



ErgoSupport® Anti-Fatigue Mat

- Intended to help reduce fatigue that healthcare personnel experience during long hours of standing in one location
- Uses a dual foam and gel technology to provide comfort under the foot
- The top surface provides a necessary friction level with damp or dry surgical booties
- The bottom surface is made of a high-traction material to provide slip resistance
- Available in four sizes
- Water-resistant
- Reusable for up to 6-12 months
- Non-sterile and not made with natural rubber latex

Part No.	Description	Color	Size	Qty/Case
82001	Anti-Fatigue Mat	Red	18" x 24"	1
82005	Anti-Fatigue Mat	Red	18" x 24"	5
83002	Anti-Fatigue Mat	Red	20" x 48"	2
84003	Anti-Fatigue Mat	Red	20" x 32"	3



Entrance Mat (TacTec)

- Intended to help remove foot borne particulate before entering the cleanroom
- Full coverage adhesive backing for adhering to the floor for use
- Seamless one-piece construction
- Numbered tabs provide easy identification and sheet separation
- Easy mobility for carts and wheeled machinery
- Polyethelene film contains a special anti-microbial agent providing long lasting protection against bacterial growth
- Non-sterile and not made with natural rubber latex

Part No.	Description	Size	Qty/Case
22-8110	Entrance Mat (TacTec)	18" x 45"	4 (30 peel off layers)



ErgoSupport® Anti-Fatigue Step Stool Mat

- Intended to help reduce fatigue that healthcare personnel experience during long hours of standing in one location
- Made with polyurethane foam to provide comfort under the foot
- Matte textured surface provides optimal friction to help prevent slips and falls
- Features a low-profile, durable edges, and slip-resistant bottom
- Sides are beveled to prevent tripping
- Small size fits on operating room step stool
- Fluid and dirt resistant with anti-microbial properties
- Non-sterile and not made with natural rubber latex

Part No.	Description	Sterility	Color	Qty/Case
85005	Anti-Fatigue Mat	Non-sterile	Black	5




ErgoSupport® Anti-Fatigue Step Stool

- Stackable and interlocking cast aluminum step stool
- Protected with an anti-microbial coating
- Designed to eliminate sites for liquid pooling and debris accumulation to help reduce contamination
- High density rubber foot pegs
- Compatible with the ErgoSupport® Anti-Fatigue Step Stool Mat
- Anti-slip surface for better stability
- Non-sterile and not made with natural rubber latex

Part No.	Description	Sterility	Color	Qty/Case
ST-001	Step Stool	Non-sterile	Silver	4





WaterBug Quiet Floor Suction® Device

Product Details

- Efficiently suctions fluid waste to help reduce risks and maintain a dry operating room floor
- Its raised design and 12' of flexible tubing allows for easy repositioning
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Tubing	Qty
90010	WaterBug Quiet Floor Suction® Device	Yes	10/Case
90020	WaterBug Quiet Floor Suction® Device	Yes	20/Case
90030	WaterBug Quiet Floor Suction® Device	Yes	30/Case



WaterBoom® Floor Suction Strip

Product Details

- Helps hold back, contain, and quietly suction an unlimited amount of waste fluid
- Provides strong suction, while conforming to your O.R. floor
- Low profile helps you step over it without suction loss
- 12' tubing attaches from the side to help reduce tripping hazards
- This unobtrusive device does not require repositioning and is a good choice for those who prefer nothing at their feet
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Length	Tubing	Qty
94010	WaterBoom® Floor Suction Strip	32"	Yes	10/Case
94030	WaterBoom® Floor Suction Strip	32"	Yes	30/Case



AquaNot® Floor Suction Device

Product Details

- Has similar characteristics as the WaterBug Quiet Floor Suction® Device, but has a lower profile
- Efficiently suctions fluid waste to help reduce risks and maintains a dry operating room floor
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Tubing	Qty
97020	AquaNot® Floor Suction Device	Yes	20/Case
97060	AquaNot® Floor Suction Device	Yes	60/Case
97050-NT	AquaNot® Floor Suction Device	No	50/Case



ViraSorb® Solidifier

Product Details

- Evenly gels canister contents in seconds without stirring or agitating
- Non-toxic, Glutaraldehyde-free, and can be safely incinerated
- It does not revert back to a fluid state, helping to eliminate spillage and leakage, qualifying it as Solid Regulated Medical Waste
- Pour nozzles are included with multi-use bottles
- Non-sterile and not made with natural rubber latex

Part No.	Description	Solidifies	Qty
1000	ViraSorb® Solidifier – Single-Use Bottle	1,200cc	96/Case
1001	ViraSorb® Solidifier – Single-Use Bottle	1,500cc	96/Case
1002	ViraSorb® Solidifier – Single-Use Bottle	2,000cc	36/Case
1003	ViraSorb® Solidifier – Single-Use Bottle	3,000cc	36/Case
1005	ViraSorb® Solidifier – Single-Use Pouch	500cc	200/Case
1010	ViraSorb® Solidifier – Multi-Use Bottle	10,000cc	12/Case
1020	ViraSorb® Solidifier – Multi-Use Bottle	20,000cc	12/Case



- | Part No. | Description                       | Qty/Box |
|----------|-----------------------------------|---------|
| AB100    | AspenBasin® Fluid Waste Container | 10      |

[illegible]

### Notes:

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

### Notes:

[illegible]





Products By Original Manufacturer

Product #	Product Description	Aspen Part#	Aspen Description	QTY
ALLEN®				
A-90019	FlexiForm® Beach Chair Mask	711003	Indy Beach Chair Face Mask	10 per Box
A-90023	Universal™ Head Restraint	711004	Universal Head Restraint	10 per Box
A-90021	Arm Positioner Restraint Disposable	711400	Norton Arm Kit	10 per Box
A-93001	Traction Boot Disposables	712900	Lindy Traction Boot Liner	6 per Box
R-S 4133-V	Leg Holder Disposable Pads	713600	Leg Comfort Pad	12 per Box
ARTHREX®				
AR-1000	Transosseous Suture Passer	710000	Hoffee® Blue Lasso Suture Retriever	5 per Box
AR-1560-01	PivotPost® Pad	712500	Parent Post Knee Pad	5 per Box
AR-1606	STAR® Sleeve	711500	Dr. Don’s® Arm Canoe	5 per Box
AR-1627-06	Universal Head Restraint	711004	Universal Head Restraint	10 per Box
AR-1635	Lateral Arm Traction Sleeve	711504	Dr. Don’s® Arm Suspension Kit – Northern	10 per Box
AR-1644	Trimano® Beach Chair Kit	711200	Berkette Beach Chair Limb Positioner Kit	6 per Box
AR-1712	Ankle Distraction Strap	713400	Smoot Ankle Distraction Strap	10 per Box
AR-6529-11	HDS Disposable Hip Kit	712100	Lewiston Positioner Kit	1 per Box
DEROYAL®				
32-525	Shoulder Suspension Kit	711502	Dr. Don’s® Arm Suspension Kit	10 per Box
IMP®				
1092P	Universal Lateral Positioner® Pads	712600	Lateral Positioner Pad	15 per Box
294P	De Mayo Knee Positioner® Pad	713100	Hermanator® Knee Positioning Pad	10 per Box
803-GP-10	De Mayo Knee Positioner® Pad	713300	XL Hermanator® Pink Knee Positioning Pad	10 per Box
803-P-10	De Mayo Knee Positioner® Pad	713100	Hermanator® Knee Positioning Pad	10 per Box
907-GP-10	Pressure Protector Pad® Set	713200	TLC Distraction Pad Set	10 per Box
INSTRAMED®				
IN-9300CLF	Upper Extremity Drape	711502	Dr. Don’s® Arm Suspension Kit	10 per Box
MCCONNELL				
12-401	Arm Support	711300	Jessie Sterile Arm Suspension Kit	10 per Box
MIZUHO OSI®				
6855-13	Hana® Perineal Post Cover	712400	Perineal Post Cover	12 per Box
5823-10	Mizuho OSI® Laminectomy Arm Cradle	714200	Sedona Laminectomy Arm Cradle	12 per Box



Products By Original Manufacturer

Product #	Product Description	Aspen Part#	Aspen Description	QTY
SCHUREMED®				
#508-0103	Disposable Schure Grip	713600	Leg Comfort Pad	12 per Box
#800-0039	Disposable Head & Chin Strap	711005	Beach Chair Head and Chin Strap	20 per Box
#800-0164	Disposable Face Mask	711000	Bethesda Mask	10 per Box
SKYTRON®				
3-040-22	Beach Chair Head & Chin Straps	711005	Beach Chair Head and Chin Strap	20 per Box
SMITH & NEPHEW®				
3757	Disposable Leg Holder pads	713600	Leg Comfort Pad	12 per Box
7111-1579	HEWSON® Suture Retriever	710000	Hoffee® Blue Lasso Suture Retriever	5 per Box
014407	GUHL™ Ankle Distraction Strap	713400	Smoot Ankle Distraction Strap	10 per Box
7210558	T-MAX™ Replacement Strap	711002	Beach Chair Head Rest Replacement Straps	1 pair per Box
7210559	T-MAX™ Disposable Face Mask	711001	Erin Mask	10 per Box
7210573	SPIDER® Shoulder Stabilization Kit	711100	Marco Shoulder Stabilization Kit	12 per Box
72200195	Shoulder Suspension Kit	711503	Dr Don’s® Arm Suspension Kit – Montana	10 per Box
72200634	Perineal Post Pad	712200	Carson 8" Perineal Post Pad	1 per Box
72200635	Boot Insert Pad	712300	Patient Care Boot Insert Pad	1 Pair per Box
72201426	SPIDER® Boot Liner	712300	Patient Care Boot Insert Pad	1 Pair per Box
72201813	FERKEL® Thigh System Pad	713500	Thigh Holder Pad	6 per Box
STERIS®				
BF039	Head Restraint Straps	711005	Beach Chair Head and Chin Strap	20 per Box
BF049	Traction Boot Liner	712900	Lindy Traction Boot Liner	6 per Box
STRYKER®				
OZ 7511975	Arthroscopy Leg Holder Pad	713600	Leg Comfort Pad	12 per Box
OZ 7518375	Chin Strap	711005	Beach Chair Head and Chin Straps	20 per Box
OZ 7547575	Disposable Mask	711000	Bethesda Mask	10 per Box
ZIMMER BIOMET®				
1320-015	Alvarado Disposable Liner	713100	Hermanator® Knee Positioning Pad	10 per Box
905026	Shoulder Traction Kit	711503	Dr. Don’s® Arm Suspension Kit-Montana	10 per Box

Allen, C-Flex, Flexiform, and Universal are either registered trademarks or trademarks of Allen Medical Systems, Inc. ("Allen Medical"). Arthrex, PivotPost, and S.T.a.R. are either registered trademarks or trademarks of Arthrex, Inc. ("Arthrex"). Trimano is a registered trademark of MAQUET GmbH ("MAQUET"). DeRoyal is a registered trademark of DeRoyal Industries, Inc. ("DeRoyal"). IMP, Innovative Medical Products, Inc., McGuire, Universal Lateral Positioner, Knee Positioning Triangles, De Mayo Knee Positioner, De Mayo Universal Distractor, and Pressure Protector Pad, are either registered trademarks or trademarks of Innovative Medical Products, Inc. ("IMP"). InstraMed is a trademark of InstraMed, Inc. ("InstraMed"). Affiniti, Affiniti 50, Affiniti 70, ClearVue, ClearVue C9-4v, ClearVue 350, ClearVue 550, ClearVue 650, ClearVue 850, HD7 XE, HD9, HD11, HD115, and Sparq are either registered trademarks or trademarks of Koninklijke Philips N.V. ("Philips"). Mizuho OSI, Hana, ProFX, and SSXT are either registered trademarks or trademarks of Mizuho OSI, Inc. ("Mizuho OSI"). SchureMed is a trademark of Schuerch Corporation ("Schuerch"). ACUSON, ACUSON X500, ACUSON X100, ACUSON X150, ACUSON X300, ACUSON NX, ACUSON NX2, ACUSON NX3, ACUSON NX3 Elite, ACUSON P500, ACUSON S2000, Sonoline G20, Sonoline G40, Sonoline G50S, and Sonoline Antares are either registered trademarks or trademarks of Siemens Medical Solutions USA, Inc. ("Siemens"). Skytron is registered trademark of Skytron, LLC. ("Skytron"). Smith & Nephew, GUHL, T-MAX, SPIDER, SPIDER2, FERKEL, and SWANSON SUTURE PASSER are either registered trademarks or trademarks of Smith & Nephew PLC ("Smith & Nephew"). HEWSON is a registered trademark of Smith & Nephew, Inc. ("Smith & Nephew"). Steris is a registered trademark of Steris Inc. ("Steris"). Stryker is a trademark of Stryker Corporation ("Stryker"). Velcro is a registered trademark of Velcro IP Holdings LLC. Zimmer is a registered trademark of Zimmer, Inc. ("Zimmer").

Aspen Surgical Products, Inc. is not affiliated with, endorsed by, or sponsored by Allen Medical, Arthrex, DeRoyal, IMP, InstraMed, Mizuho OSI, Philips, Schuerch, Siemens, Skytron, Smith & Nephew, Steris, Stryker, or Zimmer.

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2023 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1045824 rev 4 14-MAR-2023 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

USA 616-698-7100  
Canada 888-364-7004  
[www.aspensurgical.com](http://www.aspensurgical.com)

