

Instructions for use

Web Guards

Intended Use:

The Web Guards are intended to protect surgical instruments through handling and sterilization.

Intended User:

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Instructions for Use:

- Place guard over clean instrument.
- Guards may be left on during handling, storage, and sterilization.
- Sterilize with low heat, plasma sterilization process only, such as the Sterrad® Sterilization System.

Disposal:

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Sterrad is a TRADEMARK of Johnson & Johnson Medical Inc.



Aspen Surgical Products, Inc.

6945 Southbelt Dr. SE, Caledonia, MI 49316

phone +1 616.698.7100 toll-free +1 888.364.7004

aspensurgical.com

Made in USA



Rx ONLY